

Shiitake Mushroom Extract – Technical Data Sheet (40% Polysaccharides)

Lentinus Edodes Extract
Lentinus edodes (Berk.) Sing
40% Polysaccharides (UV)
Light brown to brown powder
Characteristic mushroom-like, earthy odor
~10-15 g/L at 50-60 °C with stirring
Practically insoluble
 Disperse in warm water (~50–60 °C) Stir for 5–10 minutes Filter if necessary to remove particulate matter
Skin conditioning, soothing, antioxidant
0.5–2% (in water phase)
Store in a cool, dry, dark place in tightly sealed packaging
24 months if properly stored and unopened
1000 mg/kg bw/day (oral, rat, polysaccharide-rich extract)
>5000 mg/kg bw (oral, mouse)
No genotoxic effects observed in standard toxicological evaluations